

# BLACKline

ULTRAmicro + ULTRAmicro flex

MEGAmicro + MEGAmicro flex

The MICRO mirrors are excellent for microsurgery, featuring a thickness of just 0.6 mm.

MICROflex comes with a flexible handle for individually adjusted mirror positioning.

The front surface coating results in an exceptionally sharp reflection, which is essential in microsurgery.

The ultra-bright HD-quality mirror offers a brilliantly color-true image. This combination makes the smallest details visible and allows for optimized work.



The BLACKline-finish prevents disturbing reflections for the eyes.



# BLACKline



ULTRAmicro  
MEGAmicro



ULTRAmicro flex  
MEGAmicro flex

All MICRO mirrors have an SS-thread and are compatible with all handles that have a SS-thread.

BLACKline	1	1	2	2
each	ULTRAmicro	MEGAmicro	ULTRAmicro flex	MEGAmicro flex
Ø 3 mm	B6181	B6200	B5181	B5200
Ø 5 mm	B6182	B6220	B5182	B5220
oval 3x6 mm	B6183	B6232	B5183	B5232
drop-shaped 3x6 mm	B6184	B6240	B5184	B5240

BLACKline ERGOtouch handles/ each	SS-thread
	B442