

RELAX

RELAX – created with a perfect blend of design and ergonomics.

Exact sharpness

by Front Surface mirror coating



Highest reflection

with crystal clear brightness, especially with ULTRA FS

Comfortably light weight

for relaxed dental operations

Comfortable to hold

due to the ergonomic design with grip recess

Keeping quality

colour stable, resistant against acids as well as against plaque-indicators

Hygienic

due to the plane transition from mirror surface to mirror casing without any gap

No distracting reflections

due to the satin-finished surfaces

The instrument is autoclavable, acid-resistant and suitable for thermal disinfectors.



RELAX



The split mirror shows the crystal clear ultra-bright reflection of the ULTRA FS compared to Rhodium.

Colours	Size	Ø in mm	Content	FS Rhodium REF	FS Ultra REF
all colours	4	22	10 Asstd.	7100	7300
white	4	22	10 Pack	7101	7301
pastel yellow	4	22	10 Pack	7102	7302
pastel blue	4	22	10 Pack	7103	7303
pastel mint	4	22	10 Pack	7104	7304
pastel lavender	4	22	10 Pack	7105	7305
anthracite	4	22	10 Pack	7106	7306
pink	4	22	10 Pack	7107	7307
all colours	5	24	10 Asstd.	7110	7310
white	5	24	10 Pack	7111	7311
pastel yellow	5	24	10 Pack	7112	7312
pastel blue	5	24	10 Pack	7113	7313
pastel mint	5	24	10 Pack	7114	7314
pastel lavender	5	24	10 Pack	7115	7315
anthracite	5	24	10 Pack	7116	7316
pink	5	24	10 Pack	7117	7317