




MEGA FS Rhodium

Form follows function –
it's ergonomically designed

More comfortable cheek retraction
due to the flat, broad-shaped stem.

Silk-matt stainless steel
prevents unwanted light reflections and gives
the MEGA a high-quality finish.

Exact sharp, high-contrast mirror image
down to the smallest detail.



Easy cleaning
due to the anti-stick effect of the special
HAHNENKRATT front surface coating.



MEGA FS Rhodium



MEGA FS Rhodium					SS-thread REF	CS-thread* REF
1	MEGA	4 plane	Ø 22 mm	box of 12	6060	6070
	MEGA	5 plane	Ø 24 mm	box of 12	6061	6071
2	MEGAduo, double sided	4 plane	Ø 22 mm	box of 6	6100	6170
	MEGAduo, double sided	5 plane	Ø 24 mm	box of 6	6101	6171
3	MEGAretract	4 plane	Ø 22 mm	box of 6	6660	6665
	MEGAretract	5 plane	Ø 24 mm	box of 6	6661	6666
4	MEGAretract, bifurcated stem	4 plane	Ø 22 mm	box of 6	6680	6685
	MEGAretract, bifurcated stem	5 plane	Ø 24 mm	box of 6	6681	6686
5	MEGAstretch		18x34 mm	each	6860	
	MEGAstretch		22x38 mm	each	6861	
6	MEGA Parallelometer		22x37 mm	each	6900	
	MEGA Wangenabhalter	without mirror	Ø 22 mm	each	6070 W	
	MEGAabrasiv <small>Saphire-hard glass surface suitable for use with aluminium oxide airbrasive</small>	4 plane	Ø 21 mm	each	6080	

***) Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.**