


MEGA FS

The alternative

The MEGA FS TiCr is the alternative to the MEGA FS Rhodium.

Highly reflective layers of titanium and chromium create a front surface mirror that reflects in a cooler hue compared to FS Rhodium.

Precisely sharp, high-contrast reflection down to the finest detail.



Easy cleaning

due to the anti-stick effect of the special HAHNENKRATT front surface coating.

Unlimited resistance

against disinfection and sterilisation, even against acid containing refresh liquids.



MEGA FS



MEGA FS / box of 12		SS-thread	CS-thread ^{*)}
size	Ø in mm	plane	plane
		REF	
4	22	6560	6565
5	24	6561	6566

^{*)} Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.