

TOPvision FS Rhodium

Front Surface

Images of exact sharpness and contrast
possible by one single reflection on the front surface.



Bright reflection in natural colours
thanks to an optimal reflection spectrum.


High sight comfort
good for fatigue-proof working

Easy cleaning
due to the anti-stick effect of the special
HAHNENKRATT front surface coating

Unlimited chemical resistance
even against acid containing refresh liquids.



TOPvision FS Rhodium



TOPvision FS Rhodium			plane	magnifying	plane	magnifying
Box of 1 size	Box of 12 size	Ø in mm	Simple Stem (SS)		Cone Socket (CS)*	
			REF			
06		6	719X06			
08		8	719X08			
10		10	719X10			
12		12	719X12			
14		14	719X14			
16		16	719X16			
	2	18	720X2		730X2	
	3	20	721X3		731X3	
	4	22	721X4	726X4	731X4	736X4
	5	24	721X5	726X5	731X5	736X5
	6	26	720X6		730X6	
	7	28	720X7		730X7	
	8	30	720X8		730X8	

***) Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.**