

MICRO


Stainless Steel, satin finished

BLACKline coating

The micro mirrors are excellent for microsurgery, featuring a thickness of just 0.6 mm.

MICROflex comes with a flexible handle for individually adjusted mirror positioning.

The front surface coating results in an exceptionally sharp reflection, which is essential in microsurgery.



The ultra-bright HD-quality mirror offers a brilliantly color-true image. This combination makes the smallest details visible and allows for optimized work.

The satin-finished surface of the stainless steel instruments prevents distracting reflections for the eye. The same applies to the deep black, matte BLACKLINE coating.



MICRO



each	ULTRAmicro	MEGAmicro	ULTRAmicro flex	MEGAmicro flex
Ø 3 mm	6181	6200	5181	5200
Ø 5 mm	6182	6220	5182	5220
oval 3x6 mm	6183	6232	5183	5232
drop-shaped 3x6 mm	6184	6240	5184	5240



each	ULTRAmicro	MEGAmicro	ULTRAmicro flex	MEGAmicro flex
Ø 3 mm	B6181	B6200	B5181	B5200
Ø 5 mm	B6182	B6220	B5182	B5220
oval 3x6 mm	B6183	B6232	B5183	B5232
drop-shaped 3x6 mm	B6184	B6240	B5184	B5240

Info: All MICRO mirrors have an SS thread and are compatible with all handles that have a SS thread.