

ERGOform 134°C

Ergonomical design

The ERGOform Handles are lightweight and wellbalanced due to the slim shaped end of the Handle. The wider diameter in the front and especially the grip recesses provide relaxed grip and thus relaxed working.

Pleasant haptics

The satin-finished surface and the smoothly rounded edges provide pleasant haptics and handling.

Keeping quality

colour stable, resistant against acids as well as against plaque-indicators

Autoclavable up to 134°C





ERGOform 134°C, Handles for Mouth Mirrors, Length 135 mm			SS-thread REF	CS-thread ^{*)} REF
	white	each	540	560
	pastel yellow	each	541	561
	pastel blue	each	542	562
	pastel mint	each	543	563
	pink	each	544	564
	pastel lavender	each	545	565
	anthracite	each	546	566
	ocean blue	each	547	567
	turquoise green	each	548	568
	berry	each	549	569

***) Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.**