

# ERGOtouch

## Explorers

The handles are made out of satin-finished stainless steel are hollow and thus ensure low-fatigue work. The grip recess and the design with maximum grip allows safe holding while the satin-finished surface eliminates disturbing glares.

The coloured handle ends are made out of glassfiberreinforced plastic material that is successfully used in practice since 2015. The material is autoclavable, colour- and shape-stable as well as resistant against acids and plaque-indicators.

The explorer part is made out of surgical steel according to EN ISO 7492 and thus ensures a high bending strength and flexibility.

The handles are also available for mouth mirrors.

**The Explorers are permanently fixed.**

**Autoclavable up to 134°C**





ERGOtouch duo Explorers (permanently fixed)



ERGOtouch Explorers, one-sided (permanently fixed)



1 Pack.	Handles for Mouth Mirrors (SS thread)	Form A 010 Fig. 23	Form C 030 Fig. 8	Form D 040 Fig. 9	Form F 060 Fig. 6	Form G 070 Fig. 17
<b>REF</b>						
<b>ERGOtouch duo</b>						
stainless steel	667	667A+...	667C+...	667D+...	667F+...	667G+...
ERGOtouch duo can be combined with every probe tip, e.g. REF 667A+C / 667A+D						
<b>ERGOtouch</b>						
stainless steel	442	600A	600C	600D	600F	600G
white	620	620A	620C	620D	620F	620G
pastel yellow	621	621A	621C	621D	621F	621G
pastel blue	622	622A	622C	622D	622F	622G
pastel mint	623	623A	623C	623D	623F	623G
pink	624	624A	624C	624D	624F	624G
pastel lavender	625	625A	625C	625D	625F	625G
anthracite	626	626A	626C	626D	626F	626G