

ERGOform

Explorers

The handles are made out of glassfiber-reinforced plastic material. This is successfully used in practice since 2015 with our RELAX Mouth Mirrors.

The material is autoclavable, colour- and shape-stable as well as resistant against acids and plaque-indicators.

Because of the lightweight material low-fatigue work is ensured and the wider diameter in the front and especially the grip recess provide a relaxed grip.

The silk-matt surface is low-reflective and has a pleasant feel.

The explorer parts are made out of surgical steel according to EN ISO 7492 and thus ensures a high bending strength and flexibility.

The handles are also available for mouth mirrors.

The Handles are available in the same colours as the RELAX Mouth Mirrors.

The Explorers are permanently fixed.

Autoclavable up to 134°C





ERGOform Explorers, one-sided (permanently fixed)

ERGOform Handles and Explorers / each						
1 Pack.	 Handles for Mouth Mirrors (SS thread)	 Form A 010 Fig. 23	 Form C 030 Fig. 8	 Form D 040 Fig. 9	 Form F 060 Fig. 6	 Form G 070 Fig. 17
REF						
white	540	540A	540C	540D	540F	540G
pastel yellow	541	541A	541C	541D	541F	541G
pastel blue	542	542A	542C	542D	542F	542G
pastel mint	543	543A	543C	543D	543F	543G
pink	544	544A	544C	544D	544F	544G
pastel lavender	545	545A	545C	545D	545F	545G
anthrazit	546	546A	546C	546D	546F	546G
ERGOform, stainless steel, hollow						
	492					