

# Matrixes

## Stainless Steel

Micro thin:  
0,03 / 0,04 / 0,05 mm

### Outstanding tensile strength

HAHNENKRATT matrixes and matrix bands are produced in various shapes. The stainless steel offers outstanding tensile strength.

### Practical dispenser

for matrix bands on a roll. The top opening of the dispenser allows for easy sliding out of the matrix band, and if necessary, the band can also be easily pushed back into the dispenser.



Nyström 0,04 mm



Drawings on a scale of 1:1

# Matrixes

## Ivory 0,04 mm



## Tofflemire 0,04 mm



Drawings on a scale of 1:1

# Matrixes



Matrixes 0,04 mm	REF						
box of 12	40X1	40X2	40X3	40X4	40X5	40X6	40X7
	40X11	40X12	40X13	40X14	40X15	40X16	
	40X101	40X102	40X103	40X113	40X150	40X151	
box of 100	41X1	41X2	41X3	41X4	41X5	41X6	41X7
	41X11	41X12	41X13	41X14	41X15	41X16	
	41X101	41X102	41X103	41X113	41X150	41X151	

Matrix bands in cardboard dispenser										
roll of 1m thickness in mm			width	roll of 3 m thickness in mm			width	roll of 10m thickness in mm		
0,03	0,04	0,05		0,03	0,04	0,05		0,03	0,04	0,05
10134	10144	10154	4 mm	10334	10344	10354	4 mm	10034	10044	10054
10135	10145	10155	5 mm	10335	10345	10355	5 mm	10035	10045	10055
10136	10146	10156	6 mm	10336	10346	10356	6 mm	10036	10046	10056
10137	10147	10157	7 mm	10337	10347	10357	7 mm	10037	10047	10057