

# ERGOform

Stainless Steel, satin finished

PEEK

## Ergonomic design

These handles are lightweight and balanced in the hand due to their slim end. The larger diameter and grip recesses allow for a comfortable grip and relaxed working.

## Pleasant haptic

The micro-finely textured, satin-matt surface and gently rounded edges provide a very pleasant feel.

## Satin-finished surfaces

These prevent unwanted light reflections and give the handle a refined finish.

## The ERGOform handle, made of satin-finished stainless steel

is hollow and feels well-balanced and lightweight in the hand.

## The ERGOform handles made of PEEK

PEEK is dimensionally and color-stable up to 210°C and can be sterilized both in an autoclave and with hot air.





Handles for Mouth Mirrors, Length 135 mm			SS-thread	CS-thread <sup>*)</sup>
ERGOform, PEEK			REF	
	ocean blue	each	580	590
	creme white	each	581	591
	turquoise	each	582	592
	violet	each	584	594
	anthracite	each	585	595
ERGOform, Stainless Steel, hollow				
		each	492	

<sup>\*)</sup> Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.