


HAHNENKRATT PSI Parodontometer

Präzise messen. Dauerhaft einsetzen. Maximaler Patientenkomfort.

REV 00, 2026-09 / AG

- **Optimierte Lesbarkeit der Skalierung**
Tiefschwarze DLC-Beschichtung für maximalen Kontrast, die im Vergleich zu gelaserten Markierungen nicht verblasst
- **Dauerhaft zuverlässige Anwendung im Praxisalltag**
da absolut resistent gegenüber allen Arten von Sterilisation und Desinfektion, bis hin zu säurehaltigen ReFresh Produkten
- **Einfache Aufbereitung**
durch die fugenfreien Übergänge im Messbereich
- **Gewebeschonende und atraumatische Sondierung**
durch feine, schlank geformte Sondenteile mit sanft abgerundeter Spitze und fugenfreien Übergängen



Unsere Parodontometer gibt es in elf verschiedenen Skalierungen sowie mit unterschiedlichen Griffen. Alle Griffe sind auch für Mundspiegel erhältlich.

Die Parodontometer mit Gewinde (Arbeitsteile) passen mit dem genormten Gewinde * auf alle Griffe für Mundspiegel.

E. HAHNENKRATT GmbH

Sibylle Schepperheyn
- PRRC -

Anlagen:
Katalogseite

* Simple Stem (SS)



www.hahnenkratt.com

E. HAHNENKRATT GmbH
Dentale Medizintechnik
Benzstraße 19
DE-75203 Königsbach-Stein


Fon +49 7232 3029-0
info@hahnenkratt.com
www.hahnenkratt.com

Geschäftsführer:
Peter Schepperheyn
Dipl. Wirt.-Ing. (FH)
Hendrik Schepperheyn
Dipl. Wirt.-Ing.
Amtsgericht Mannheim
HRB 501452
DE 811 122 146

HAHNENKRATT PSI Periodontometer
Measure with precision. Sustained performance.
Designed for maximum patient comfort.

REV 00, 2026-09 / AG

- **Optimized scale readability**
 Deep-black DLC coating ensures maximum contrast and does not fade over the time, contrary to laser-engraved markings
- **Sustained reliable performance in daily clinical practice**
 completely resistant against all common sterilization and disinfection methods, including acidic ReFresh Products
- **Easy Reprocessing**
 due to seamless transitions in the measurement area
- **Gingiva tissue-preserving and atraumatic probing**
 due to slim, finely designed probe tips with softly rounded ends and seamless transitions



Our Periodontometers are available in eleven different scales and with various Handles. All Handles are also available for Mouth Mirrors.

The Periodontometers with thread (working parts) fits with their standardized thread * on all Handles for Mouth Mirrors.



www.hahnenkratt.com

E. HAHNENKRATT GmbH

Sibylle Schepperheyn
 - PRRC -

Enclosure:
 Catalog Page

* Simple Stem (SS)

E. HAHNENKRATT GmbH
 Dentale Medizintechnik
 Benzstraße 19
 DE-75203 Königsbach-Stein

Fon +49 7232 3029-0
info@hahnenkratt.com
www.hahnenkratt.com

Geschäftsführer:
 Peter Schepperheyn
 Dipl. Wirt.-Ing. (FH)
 Hendrik Schepperheyn
 Dipl. Wirt.-Ing.
 Amtsgericht Mannheim
 HRB 501452
 DE 811 122 146