

Cytec

ROOT POST SYSTEM

Cytec is designed for roots with rather cylindrical shape.

Superior Benefits:

High bending strength

HT-Glassfiber: 1.678 MPa

Determined accordingly to EN ISO 178

High fatigue and fracture resistance

Secure adhesion due to micro-retentive surface, without the use of a silane

Homogeneity due to dentin-like elasticity


Effectively and sharply cutting burs

High radiopacity

Scientifically and practically proven since 1995




For further information please see our brochure or visit our website www.hahnenkratt.com.





REF 43 600





REF 43 610



REF 10 001 (empty)

The Standard Sets and Test Sets are shipped in a SHIPPING BOX.
 Posts and Drills each in blister packs.
 Please see the price list for the contents.

| Cytec Root Post System | | | universal | 1,0 mm | 1,2 mm | 1,4 mm | 1,6 mm | 1,8 mm | 2,1 mm |
|--|-----------------------------------|--------------|------------|----------|--------|--------|----------|--------|--------|
| | | Coding | – | red | white | yellow | green | blue | black |
| | | packing unit | REF | | | | | | |
|  | Preshaping Drill with guiding tip | 1 | 42 010 | | | | | | |
|  | Preshaping Drill with cutting tip | 1 | 43 000 | | | | | | |
|  | Calibration Drill | 1 | | 4300 D10 | 43 001 | 43 002 | 4300 D16 | 43 003 | 43 004 |
| Cytec HT-Glasfiber | | Standard Set | assortment | 43 600 | | | | | |
| | | Test Set | assortment | 43 610 | | | | | |
|  | HT-Glassfiber | 10 | | 4360 D10 | 43 601 | 43 602 | 4360 D16 | 43 603 | 43 604 |
| System Box, empty | | 1 | 10 001 | | | | | | |