

Photography Mirrors

Rhodium Front Surface

ULTRA Front Surface

The precisely sharp reflection is high-contrast down to the finest detail. Combined with brilliant color accuracy and pleasantly bright reflection, you get intraoral images of great clarity.



The capture of the mirror image reveals new perspectives.



Rhodium FS
Canon 7D and Metz flash:
Aperture 25 – Shutter speed
125 – ISO 100



ULTRA FS
The unique construction of the reflective layers – consisting of 22 layer packages – provides an exceptionally bright mirror image through optimized light refraction.



Photography Mirrors



REF 6660 KO / empty



REF 6530 X ... / 6730 X ...



REF 6540 X ... / 6740 X ...



REF 6550 X ... / 6750 X ...

All mirror panels have a thickness of 3 mm.

Photography Mirrors / each		double-sided without	single-sided without handle		handle bent		handle V-shaped		handle straight		
Form	mm	REF Rhodium FS	REF Rhodium FS	REF ULTRA FS	REF Rhodium FS	REF ULTRA FS	REF Rhodium FS	REF ULTRA FS	REF Rhodium FS	REF ULTRA FS	
1		58x76	6300X1	6500X1	6700X1	6530X1	6730X1			6550X1	6750X1
		66x95	6300X5	6500X5	6700X5	6530X5	6730X5			6550X5	6750X5
2		74x106	6300X2	6500X2	6700X2	6530X2	6730X2			6550X2	6750X2
3		42x128	6300X3	6500X3	6700X3	6530X3	6730X3	6540X3	6740X3		
4		36x120		6500X3KF0	6700X3KF0	6530X3KF0	6730X3KF0	6540X3KF0	6740X3KF0		
		55/40x157	6300X4								
		53/31x142	6300X6								
		86/77x133	6300X7								
		71/62x133	6300X8								
		47x153	6300X9								
	Suitcase for Photography Mirrors / empty			REF							
For 3 Photography Mirrors with handle, empty Outside: 285x240x100 mm			6660 KO								

