

ERGO

ERGOform

HANDLES IN ERGONOMICAL DESIGN

PEEK + Stainless Steel, satin finished



ERGOform

Temperature resistant up to 210°C

The used material called **PEEK** is a plastic material which is also used for implants. **PEEK** is colour-proof and resistant to deformation up to 210°C. It can be sterilised in the autoclave, as well as with hot-air.

Ergonomical design

The ERGOform Handles are lightweight and well-balanced due to the slim shaped end of the Handle. The wider diameter in the front and especially the grip recess provide relaxed grip and thus relaxed working.

Pleasant haptics

The satin-finished surface and the smoothly rounded edges provide pleasant haptics and handling.

Satin-finished surface

eliminates disturbing glares and gives the instrument an exclusive finish.



Handles for Mouth Mirrors, Length 135 mm			SS-thread REF	CS-thread REF
ERGOform, made of PEEK Plastic 	grey-blue	each	580	590
	white	each	581	591
	turquoise-green	each	582	592
	pastel-blue	each	583	593
	aubergine	each	584	594
	anthracite	each	585	595
ERGOform, made of Stainless Steel, hollow 		each	492	