

Handles

FOR MOUTH MIRRORS

Pleasant haptics

The high quality is characterized by the careful finish. The surfaces are fine polished and all edges are smoothly rounded. Thus, the instruments provide a pleasant handling.

The hollow Handles are well-balanced and light-weight. A wider diameter of some handles and especially the grip recess provide relaxed grip and therefore relaxed working.












Satin-finished surface

eliminates disturbing glares and gives the instruments an exclusive finish

Outstanding stability against corrosion

According to our experience of over 60 years, the stainless steel is stable against corrosion. It complies with EN ISO 7153-1 (Surgical Instruments – Metallic materials).



Handles for Mouth Mirrors, Length 130 mm		SS-thread REF	CS-thread REF
Stainless Steel polished, octagonal, SE-form 	each	413	415
Stainless Steel polished, octagonal, serrated 	each	420	421
Stainless Steel polished, round, serrated 	each	426	427
Stainless Steel polished, round, serrated 	each	5012	5013
Stainless Steel polished 	each	430	
Stainless Steel polished, hollow 	each	460	461
Stainless Steel satin-finished surface, hollow 	each	466	467
Stainless Steel polished, hexagonal, hollow 	each	470	471
Stainless Steel satin-finished, hexagonal, hollow 	each	476	477
Stainless Steel polished, hexag., hollow 135 mm 	each	480	481
Stainless Steel satin-finished, hexagonal, with scale 	each	474	