

TOPvision FS

Front Surface

The TOPvision FS is the alternative to the TOPvision FS Rhodium.

Highly reflective layers of titanium and chromium create a front surface mirror that reflects in a cooler hue compared to FS Rhodium.

Precisely sharp, high-contrast reflection down to the finest detail.

only one reflection
▶ exact sharp image



Front Surface Mirror

High sight comfort

gentle to the eye and good for fatigue-proof working.

Easy cleaning


due to the anti-stick effect of the special HAHNENKRATT front surface coating.

Unlimited resistance

against disinfection and sterilisation, even against acid containing refresh liquids.



TOPvision FS



TOPvision FS / Box of 12		Simple Stem (SS)	Cone Socket (CS) *
size	Ø in mm	plane	plane
		REF	
4	22	633X4	636X4
5	24	633X5	636X5

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.

