

MEGA FS

Front Surface

The MEGA FS is the alternative to the MEGA FS Rhodium.

Highly reflective layers of titanium and chromium create a front surface mirror that reflects in a cooler hue compared to FS Rhodium.

Precisely sharp, high-contrast reflection down to the finest detail.

only one reflection
▶ exact sharp image



Front Surface Mirror

High sight comfort

gentle to the eye and good for fatigue-proof working.

Easy cleaning

due to the anti-stick effect of the special HAHNENKRATT front surface coating.

Unlimited resistance

against disinfection and sterilisation, even against acid containing refresh liquids.



MEGA FS



MEGA FS / box of 12		SS-thread	CS-thread [*]
size	Ø in mm	plane	plane
		REF	
4	22	6 560	6 565
5	24	6 561	6 566

^{*}Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
 CS is mainly used in the USA, Canada, South-Korea and Spain.

