

ULTRA FS

ULTRA FS ultra bright

FRONT SURFACE



ULTRA FS – the brightest and sharpest mirror worldwide.

The new unique design of the mirror layers provides superior bright images and crystal clarity previously unavailable.



The splitted mirror shows the crystal clear bright reflection of the ULTRA FS compared to Rhodium.

The ULTRA FS features are:

Highest reflection with crystal clear brightness


Superior colour accuracy, without any discoloration

Exact sharpness

Scratch resistant


Improved comfortable cheek retraction

only one reflection eliminates the so-called »ghost image«
→ exact sharp image



Front Surface Mirror

GRAPHIC REFLECTION



Determined by Fraunhofer Institute Germany



ULTRA FS		ATTENTION: Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.		SS-thread REF	CS-thread REF	
1	ULTRA FS	4 plane	box of 6	6 160	6 165	
	ULTRA FS	5 plane	box of 6	6 161	6 166	
2	ULTRAduo FS, double sided	4 plane	box of 6	6 640	6 645	
	ULTRAduo FS, double sided	5 plane	box of 6	6 641	6 646	
3	ULTRAretract FS	4 plane	box of 6	6 670	6 675	
	ULTRAretract FS	5 plane	box of 6	6 671	6 676	
4	ULTRAretract FS, bifurcated stem	4 plane	box of 6	6 690	6 695	
	ULTRAretract FS, bifurcated stem	5 plane	box of 6	6 691	6 696	
5	ULTRAretract FS, bifurcated stem, stem double 45°	4 plane	each	6 650		
6	ULTRAmicro	mini mirror	Ø 3 mm	each	6 181	
	ULTRAmicro	mini mirror	Ø 5 mm	each	6 182	
	ULTRAmicro	mini mirror	3 x 6 mm	each	6 183	