




HAHNENKRATT.COM

ERGOform

ERGOform EXPLORERS

ERGOform HANDLES FOR MOUTH MIRRORS



ERGOform Handles

The glassfiber-reinforced plastic material of the ERGOform Handles is successfully used in practice since 2015 with our RELAX Mouth Mirrors. The material is autoclavable, colour- and shape-stable as well as resistant against acids and plaque-indicators.

Because of the lightweight material low-fatigue work is ensured and the wider diameter in the front and especially the grip recess provide a relaxed grip. The Handles for Mouth Mirrors and Explorers are available in the same colours as the RELAX Mouth Mirrors. Thus the practice organisation can be simplified.

ERGOform Explorer

made out of surgical steel according to EN ISO 7492 ensures a high bending strength and flexibility.

The Explorers are permanently fixed.

Autoclavable up to 134°C



ERGOform | Handles and Explorers | each

The Explorers are permanently fixed.

	Handles for Mouth Mirrors REF	Form A 010 Fig. 23 REF	Form C 030 Fig. 8 REF	Form D 040 Fig. 9 REF	Form F 060 Fig. 6 REF	Form G 070 Fig. 17 REF
white	540	540 – A	540 – C	540 – D	540 – F	540 – G
pastel yellow	541	541 – A	541 – C	541 – D	541 – F	541 – G
pastel blue	542	542 – A	542 – C	542 – D	542 – F	542 – G
pastel mint	543	543 – A	543 – C	543 – D	543 – F	543 – G
pink	544	544 – A	544 – C	544 – D	544 – F	544 – G
pastel lavender	545	545 – A	545 – C	545 – D	545 – F	545 – G
anthracite	546	546 – A	546 – C	546 – D	546 – F	546 – G