

# BLACKline

The BLACKline instruments receive with a DLC (Diamondlike carbon) coating a high-end finish.

The matt finish eliminates disturbing reflections. This is an important aspect for working under the microscope, which is becoming increasingly important in endodontics.

The BLACKline increases working comfort while adding ergonomic elements and offers superior advantages for every purpose.

# BLACKline

ULTRAmicro + ULTRAmicro flex


MEGAmicro + MEGAmicro flex

The micro mirrors are excellent for microsurgery, featuring a thickness of just 0.6 mm.


MICROflex comes with a flexible handle for individually adjusted mirror positioning.

The front surface coating results in an exceptionally sharp reflection, which is essential in microsurgery.

The ultra-bright HD-quality mirror offers a brilliantly color-true image. This combination makes the smallest details visible and allows for optimized work.




The BLACKline-finish prevents disturbing reflections for the eyes.





ULTRAmicro  
MEGAmicro



ULTRAmicro flex  
MEGAmicro flex

INFO: All MICRO mirrors have an SS thread and are compatible with all handles that have a SS thread.

BLACKline	1	1	2	2
each	ULTRAmicro	MEGAmicro	ULTRAmicro flex	MEGAmicro flex
Ø 3 mm	B6181	B6200	B5181	B5200
Ø 5 mm	B6182	B6220	B5182	B5220
oval 3x6 mm	B6183	B6232	B5183	B5232
drop-shaped 3x6 mm	B6184	B6240	B5184	B5240

BLACKline ERGOtouch handles/ each	SS-thread
	B442

