

ERGOtouchtouch

ERGOtouch EXPLORERS
ERGOtouch HANDLES FOR MOUTH MIRRORS





ERGOtouch Handles

made out of satin-finished stainless steel are hollow and thus ensure low-fatigue work. The grip recess and the design with maximum grip allows safe holding while the satin-finished surface eliminates disturbing glares.

The coloured handle ends are made out of glassfiberreinforced plastic material that is successfully used in practice since 2015. The material is autoclavable, colour- and shape-stable as well as resistant against acids and plague-indicators.


ERGOtouch Explorer

made out of surgical steel according to EN ISO 7492 ensures a high bending strength and flexibility.

The Explorers are permanently fixed.

Autoclavable up to 134°C





442 (SS)

1

stainless steel



ERGOtouch Explorers, one-sided (permanently fixed)

600G

	Handles for Mouth Mirrors REF each	Form A 010 Fig. 23 REF each	Form C 030 Fig. 8 REF each	Form D 040 Fig. 9 REF each	Form F 060 Fig. 6 REF each	Form G 070 Fig. 17 REF each
stainless steel	667	667 A +	667 C +	667 <mark>D+</mark>	667F+	667 <mark>G+</mark>

ERGOtouch duo can be combined with every probe tip, e.g. REF 667A+C / 667A+D

Stanness steel	443 (CS)	00071	0000	0000	0001	0000				
SS: Simple Stem CS: Cone Socket										
white	620	620A	620C	620D	620F	620G				
pastel yellow	621	621A	621C	621D	621F	621G				
pastel blue	622	622A	622C	622D	622F	622G				
pastell mint	623	623A	623C	623D	623F	623G				
pink	624	624A	624C	624D	624F	624G				
pastel lavender	625	625A	625C	625D	625F	625G				
anthracite	626	626A	626C	626D	626F	626G				