

TOPvision FS

TOPvision FS Rhodium


FRONT SURFACE RHODIUM

TOPvision FS Rhodium are coated with Rhodium on the front surface. Rhodium is a precious metal as gold and as a precious metal it's chemical resistant, even against acids.

The TOPvision FS Rhodium features are:


Images of exact sharpness and contrast up to the smallest detail, without any distortion

only one reflection
→ exact sharp image



Front Surface Mirror

Two reflections
→ double image



Rear Surface Mirror


Bright and brilliant reflection in natural colours
thanks to an optimal reflection spectrum.

Highest sight comfort
gentle to the eye and good for fatigue-proof working

Easy cleaning
due to the anti-stick effect of the special
HAHNENKRATT front surface coating

Unlimited resistance
against disinfection and sterilisation, even against
acid containing refresh liquids.





TOPvision FS Rhodium mouth mirrors			REF			
Box of 1 size	Box of 12 size	Ø in mm	Simple Stem (SS)		Cone Socket (CS)	
			plane	magnifying	plane	magnifying
06		6	719X06			
08		8	719X08			
10		10	719X10			
12		12	719X12			
14		14	719X14			
16		16	719X16			
	2	18	720X2		730X2	
	3	20	721X3		731X3	
	4	22	721X4	726X4	731X4	736X4
	5	24	721X5	726X5	731X5	736X5
	6	26	720X6		730X6	
	7	28	720X7		730X7	
	8	30	720X8		730X8	