



HAHNENKRATT.COM

Economy

MOUTH MIRRORS



THE ECONOMICALLY PRICED

High image quality


The Economy mouth mirrors have a threefold pure silver coating. The reflection is bright and offers high image quality.

Best stability against corrosion

According to our experience of over 50 years, the stainless steel is stable against corrosion. It complies with ISO 7153-1 (Surgical Instruments – Metallic materials).

High resistance against disinfection and sterilisation

The silver coating is protected with a copper layer and it's additionally sealed with a special protective film.



Economy mouth mirrors		REF			
Box of 12 size	Ø in mm	Simple Stem (SS)		Cone Socket (CS)	
		plane	magnifying	plane	magnifying
4	22	770X4	771X4	775X4	776X4
5	24	770X5	771X5	775X5	776X5