

# Handles

## FOR MOUTH MIRRORS

### **Pleasant haptics**

The high quality is characterized by the careful finish. The surfaces are fine polished and all edges are smoothly rounded. Thus, the instruments provide a pleasant handling.

The hollow Handles are well-balanced and light-weight. A wider diameter of some handles and especially the grip recess provide relaxed grip and therefore relaxed working.













### **Satin-finished surface**

eliminates disturbing glares and gives the instruments an exclusive finish

### **Outstanding stability against corrosion**

According to our experience of over 50 years, the stainless steel is stable against corrosion. It complies with ISO 7153-1 (Surgical Instruments – Metallic materials).



| Handles for Mouth Mirrors, Length 130 mm  |      | SS-thread<br>REF | CS-thread<br>REF |
|---|------|------------------|------------------|
| Stainless Steel polished, octagonal, SE-form<br>   | each | 413              | 415              |
| Stainless Steel polished, octagonal, serrated<br>  | each | 420              | 421              |
| Stainless Steel polished, round, serrated<br>  | each | 426              | 427              |
| Stainless Steel polished, round, serrated<br>  | each | 5012             | 5013             |
| Stainless Steel polished<br>   | each | 430              | 431              |
| Stainless Steel polished, hollow<br>   | each | 460              | 461              |
| Stainless Steel satin-finished surface, hollow<br>   | each | 466              | 467              |
| Stainless Steel polished, hexagonal, hollow<br>  | each | 470              | 471              |
| Stainless Steel satin-finished, hexagonal, hollow<br>                                      | each | 476              | 477              |
| Stainless Steel polished, hexag., hollow 135 mm<br>  | each | 480              | 481              |
| Chrome-plated brass, hexagonal, hollow<br>   | each | 422              | 423              |
| Stainless Steel satin-finished, hexagonal, with scale<br>(for orientation measurement)<br> | each | 474              |                  |