

Contec

ROOT POST SYSTEM

Contec is designed for roots with rather conical shape. The shape is designed with two different tapers and follows the natural shape of the root canal.

Overview of your advantages:

Utmost bending strength

HT-Carbonfiber: 1.857 MPa

HT-Glassfiber: 1.678 MPa

Determined accordingly to EN ISO 178

Utmost fatigue and fracture resistance


Homogeneity due to dentin-like elasticity

Safe adhesion due to micro-retentive surface

Effectively and sharply cutting burs

High radiopacity

Scientifically and practically proved since 1995



For further information
please see our brochure
or visit our website
www.hahnenkratt.com.






REF 44 600








REF 44 610



REF 10 002 + REF 10 000

NEW

Contec Root Post System		Code	universal	1,1 mm	1,3 mm	1,5 mm	1,75 mm	2,0 mm
			–	red	white	yellow	blue	black
		package of	REF					
	Preshaping Drill with guiding tip	1	42 010					
	Preshaping Drill with cutting tip	1	43 000					
	Calibration Drill	1		44 00 D1,1	44 001	44 002	44 003	44 004
Contec Carbon HT-Carbonfiber	Set Standard	Assortment	44 400					
	Trial Kit	Assortment	44 410					
	HT-Carbonfiber	10			44 401	44 402	44 403	44 404
Contec blanco HT-Glasfiber	Set Standard	Assortment	44 600					
	Trial Kit	Assortment	44 610					
	HT-Glasfiber	10		44 60 D1,1	44 601	44 602	44 603	44 604
System Box + Organizer, empty		1	10 002 + 10 000					

Content of the assortment: Please see our price list.