




HAHNENKRATT.COM

IDS News 2017

# ERGOform

ERGOform EXPLORERS

ERGOform HANDLES FOR MOUTH MIRRORS



## ERGOform Handles

The glassfiber-reinforced plastic material of the ERGOform Handles is successfully used in practice since 2015 with our RELAX Mouth Mirrors. The material is autoclavable, colour- and shape-stable as well as resistant against acids and plaque-indicators.

Because of the lightweight material low-fatigue work is ensured and the wider diameter in the front and especially the grip recess provide a relaxed grip. The Handles for Mouth Mirrors and Explorers are available in the same colours as the RELAX Mouth Mirrors. Thus the practice organisation can be simplified.

## ERGOtouch Explorer

made out of surgical steel according to EN ISO 7492 ensures a high bending strength and flexibility.













**The Explorers are permanently fixed.**

**Autoclavable up to 134°C**



**ERGOform** | Handles and Explorers | each

**The Explorers are permanently fixed.**

							
	<b>Handles for Mouth Mirrors</b> REF	<b>Form A</b> 010 Fig. 23 REF	<b>Form C</b> 030 Fig. 8 REF	<b>Form D</b> 040 Fig. 9 REF	<b>Form F</b> 060 Fig. 6 REF	<b>Form G</b> 070 Fig. 17 REF	
	white	540	540 – A	540 – C	540 – D	540 – F	540 – G
	pastel yellow	541	541 – A	541 – C	541 – D	541 – F	541 – G
	pastel blue	542	542 – A	542 – C	542 – D	542 – F	542 – G
	pastel mint	543	543 – A	543 – C	543 – D	543 – F	543 – G
	pink	544	544 – A	544 – C	544 – D	544 – F	544 – G
	pastel lavender	545	545 – A	545 – C	545 – D	545 – F	545 – G
	anthracite	546	546 – A	546 – C	546 – D	546 – F	546 – G