





HAHNENKRATT.COM

# ERGOtouch

## ERGOtouch

ERGOtouch EXPLORERS

ERGOtouch HANDLES FOR MOUTH MIRRORS



### ERGOtouch Handles

made out of satin-finished stainless steel are hollow and thus ensure low-fatigue work. The grip recess and the design with maximum grip allows safe holding while the satin-finished surface eliminates disturbing glares.

The coloured handle ends are made out of glassfiber-reinforced plastic material that is successfully used in practice since 2015. The material is autoclavable, colour- and shape-stable as well as resistant against acids and plaque-indicators.

### ERGOtouch Explorer

made out of surgical steel according to EN ISO 7492 ensures a high bending strength and flexibility.

**The Explorers are permanently fixed.**






**Autoclavable up to 134°C**



ERGOtouch duo Explorers (permanently fixed)



ERGOtouch Explorers, one-sided (permanently fixed)

			 <p><b>Form A</b> 010 Fig. 23 REF   each</p>	 <p><b>Form C</b> 030 Fig. 8 REF   each</p>	 <p><b>Form D</b> 040 Fig. 9 REF   each</p>	 <p><b>Form F</b> 060 Fig. 6 REF   each</p>	 <p><b>Form G</b> 070 Fig. 17 REF   each</p>
	<b>Handles for Mouth Mirrors</b>						
	REF   each						
1	stainless steel	667	667 – A + ...	667 – C + ...	667 – D + ...	667 – F + ...	667 – G + ...

ERGOtouch duo can be combined with every probe tip, e.g. REF 667 - A+C / 667 - A+D

2	stainless steel	442 (SS) 443 (CS)	600 – A	600 – C	600 – D	600 – F	600 – G
		SS: Simple Stem CS: Cone Socket					

	white	620	620 – A	620 – C	620 – D	620 – F	620 – G
	pastel yellow	621	621 – A	621 – C	621 – D	621 – F	621 – G
	pastel blue	622	622 – A	622 – C	622 – D	622 – F	622 – G
	pastel mint	623	623 – A	623 – C	623 – D	623 – F	623 – G
	pink	624	624 – A	624 – C	624 – D	624 – F	624 – G
	pastel lavender	625	625 – A	625 – C	625 – D	625 – F	625 – G
	anthracite	626	626 – A	626 – C	626 – D	626 – F	626 – G